


# LAPAROSUCK™

**Sterile Disposable Laparoscopic Surgery Tubing kit**





“ Providing the best solutions in single use sterile medical devices”

*Model : 2111*

**Specifications:**

- **Sterile & Disposable**
- **Soft, Flexible kink Resistant Medical Grade Tubing**
- **Flexible Colour Coded Connector**
- **Prevents the risk of cross contamination**
- **Luer Connector for Insufflation Cannula**
- **0.027 Micron Insufflation Bacterial / Viral Filter**
- **Models Designed as per your Convenience**
- **99.999% Filter Efficiency**



Model	Design	Description
1011		Sterile Suction, Irrigation, Insufflation tubing 2.5m, Colour coded connectors & Luer connector
2111		Sterile Suction, Insufflation tubing 2.5m, Irrigation tubing 2m with Spike & flow regulator, Colour coded connectors & Luer connector
3033		Sterile Suction, Irrigation, Insufflation tubing 2.5m, Colour coded connectors, Luer connector & 0.027 micron filter
3044		Sterile Suction, Insufflation tubing 2.5m, Irrigation tubing 2m with Spike & flow regulator, Colour coded connectors, Luer connector & 0.027 micron filter
Laparosuck Plus		Sterile Suction & Irrigation Trumpet Valve with Cannulae & Tubing and 0.027 micron filter
Laparosuck-GTF		Sterile Insufflation tubing 2.5m, Colour coded connector, Luer connector for insufflation cannula, 0.027 micron filter
Laparoguard		Sterile 0.027 micron insufflation Filter Bacterial & Viral filter.



## Anaesthesia

<b>SUCTOCATH<sup>®</sup></b>	Suction Catheter
<b>SUCTOCATH - TC<sup>®</sup></b>	Suction Catheter with thumb control
<b>OXYGEN MASK</b>	
<b>OXYGEN SET</b>	Twin Bore Nasal Oxygen Set
<b>GUEDAL AIRWAYS</b>	
<b>ET TUBE</b>	

## Urology

<b>TOPDRAIN PLUS<sup>®</sup></b>	Urine Bag 2000ml with Top outlet
<b>BOTDRAIN PLUS<sup>®</sup></b>	Urine Bag with Bottom outlet
<b>LEGUBAG<sup>®</sup></b>	Leg Urine Bag
<b>MOLIUMETER<sup>®</sup></b>	Urine Bag with Meter
<b>MOLICATH<sup>®</sup></b>	Latex / Silicone Foley Balloon Catheter
<b>KINGCATH<sup>®</sup></b>	Urethral Nelaton Catheter
<b>QUEENCATH<sup>®</sup></b>	Urethral Catheter
<b>ENDOFLOW<sup>®</sup></b>	Irrigation Set for Endoscopic T.U.R
<b>ENDOFLOW PLUS<sup>®</sup></b>	Endoscopic Irrigation Set (TUR) with Outlet Drainage Connection
<b>MOLICOND<sup>®</sup></b>	Male External Catheter
<b>PAEDIATRIC URINE COLLECTING BAG</b>	

## Gastro - Enterology

<b>RYLES TUBE</b>	
<b>INFANT FEEDING TUBE</b>	
<b>NUTRIFEED DRIP SET<sup>®</sup></b>	Enteral Feeding Set
<b>MOLIEX<sup>®</sup></b>	Mucus Extractor

Dealer :

## Cardio / General Surgery

<b>CHESTOCATH<sup>®</sup></b>	Thoracic / Chest Drainage Catheter
<b>MOLIVAC<sup>®</sup></b>	Closed Wound Suction Set
<b>MOLISUCK<sup>®</sup></b>	Yankaur Suction Set
<b>MOLI ADK<sup>®</sup></b>	Abdominal Drainage Kit
<b>MOLIDRAIN<sup>®</sup></b>	Underwater Sealed Drainage System
<b>CORRUGATED DRAINAGE SHEET</b>	
<b>LAPAROSUCK<sup>™</sup></b>	Sterile Disposable Laparoscopic Surgery Tubing Kit
<b>LAPAROSUCK PLUS<sup>™</sup></b>	Sterile Suction & Irrigation Trumpet Valve with Cannulae & Tubing and 0.027 micron Filter
<b>LAPAROSUCK-GTF<sup>™</sup></b>	Sterile Insufflation Tubing with 0.027 micron Filter
<b>LAPAROGUARD<sup>™</sup></b>	Laparoscopic Insufflation 0.027 micron Filter

## Infusion / Transfusion

<b>INFUSION SET</b>	
<b>BLOOD ADMINISTRATION SET</b>	
<b>MICRO INFUSION SET</b>	
<b>PRESSURE MONITORING LINES</b>	
<b>THREE WAY STOPCOCK</b>	
<b>THREE WAY STOP COCK WITH EXTENSION</b>	
<b>EXTENSION LINES</b>	
<b>MEASURED VOLUME FLUID ADMINISTRATION SET</b>	
<b>MOLIFLOW PLUS<sup>®</sup></b>	IV Flow Control Regulator with Infusion Set
<b>MOLIGUARD PLUS<sup>®</sup></b>	0.2 Micron Fluid Filter with Infusion Set
<b>PERITONEAL DIALYSIS TRANSFUSION SET</b>	

## Miscellaneous Products

<b>MOLICLAMP<sup>®</sup></b>	Umbilical Cord Clamp
<b>NEBULIZER KIT</b>	Mask with Nebulizer Chamber
<b>SURGEON GOWN</b>	
<b>FACE MASK</b>	
<b>SURGEON CAP</b>	
<b>MOLISAFE BED</b>	Sterile Bed Protection Sheet
<b>MOLISAFE ADULT DIAPER</b>	
<b>DRAPES</b>	
<b>MEDICAL GRADE TUBE DIFFERENT SIZES</b>	
	(SILICONE / PVC / LATEX )
<b>MOLISAFE LATEX EXAMINATION GLOVES<sup>®</sup></b>	



CE ISO G.M.P

**MORRISONS LIFECARE PVT. LTD.**

No:3 Anna Street, Padikuppam, Chennai - 600 107, India.  
Ph: 91 - 44 - 26155047 / 26156047, Telefax: 91 - 44 - 26154047  
email@morrisonslifecare.com, sales@morrisonslifecare.com

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)