

FEEDING DRIP SET

TWO WAY



Product Features :

- Sharp piercing spike for easy insertion in gravity feeding bag and container.
- Provided with dual drip chamber controlled by roller clamp for feeding two different fluids simultaneously.
- Made of non-toxic medical grade PVC easy visualizing.
- Efficient roller controller ensures better flow control and is easy to operate.
- Super smooth, clear, transparent, kink-resistant tube of 220 cm long tubing for unobstructed flow of fluid.
- Fixed male connector for secure connection to all standard devices.



ORDERING INFORMATION

Ref No	Model	Qty / Box	Qty / Case
7010101	One Way	80	320
7010102	Two Way	40	

 **MORRISON'S**
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com