


INFANT FEEDING TUBE



Product Features :

- For neonatal pediatric nutritional feeding
- Soft , frosted and kink resistant tubing
- Soft and rounded tip to prevent trauma during application
- Proximal end is fitted with female luer mount for easy connection to feeding funnel
- Colour coded for size identification

ORDERING INFORMATION

Ref No	Size	Qty / Minipack	Qty / Box
7010103	5fg	50	150
7010104	6fg	50	150
7010105	7fg	50	150
7010106	8fg	50	150
7010107	9fg	50	150
7010108	10fg	50	150
7010109	12fg	50	150



 **MORRISONS**
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com