

RYLES TUBE



Product Features :

- Suitable for nasogastric insertion for nutritional purposes.
- Made of medical grade non toxic, non-irritant PVC.
- Radioopaque line allows confirmation of tube placement.
- Four lateral eyes are provided for efficient aspiration & administration.
- The tube is marked at 50, 60 & 70cm from the tip for accurate placement.
- Proximal end is provided with universal funnel shape connector.

ORDERING INFORMATION

Ref No	Size	Qty / Box	Qty / Case
7010110	6fg	300	1200
7010111	8fg	300	1200
7010112	10fg	300	1200
7010113	12fg	250	1000
7010114	14fg	200	800
7010115	16fg	200	800
7010116	18fg	150	600
7010117	20fg	150	600



 **MORRISONS**
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com