

# MOLISYTE™

## Luer Access Split Septum



- Needle Free IV Access Device.
- Transparent housing material made of polycarbonate allows visual assessment of fluid path.
- Low priming volume : 0.8 / 0.45ml of length 15 cm.
- Sterile & Available in Singleway, Biway & Triway Extension line.
- Latex Free

PRODUCT	DESCRIPTION
Molisyte-Q	Needle Free Connector
Molisyte Single-Way	Single-Way Extension Set-15cm(6IN)-0.45ml/0.8ml
Molisyte Bi-Way	Bi-Way Extension Set-15cm(6IN)-0.45ml/0.8ml
Molisyte Tri-Way	Tri-way Extension Set-15cm(6IN)-0.45ml/0.8ml

### ORDERING INFORMATION

Ref No	Model / Size	Qty/Box	Qty/Case
1010178	Single Way	200	800
1010179	Bi Way	200	800
1010180	Tri Way	200	800
1010181	Molisyte-Q	1000	4000



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)