


LIFECARE INFUSION SET NON-VENTED



- Flow rate of 20drops/ml
- Plain Latex for intermediate medication
- Tip connector for Needle connection
- Universal Non-Vented Spike
- 150cm length kink resistant medical grade PVC tube
- Provided with 21G Needle



ORDERING INFORMATION

Ref No	Qty/Box	Qty/Case
1010103	25	500

 **MORRISONS**
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com