

# MORRISONS®

## INFUSION SET - NON VENTED

### I.V Set JSM



Bulb latex for intermediate medication

#### Product Features :

- For gravity feed only.
- Sharp piercing spike for easy insertion to IV container.
- Transparent drip chamber for easy monitoring of flow rate and priming.
- Kink resistant medical grade PVC tubing.
- Smooth roller clamp for easy flow control.
- 15 micron fluid filter for filtration of foreign particles.
- Drop rate of 20 drops /ml.
- Standard tube length-150cm.
- Single use • Disposable.
- "Sterile, individually packed"

#### ORDERING INFORMATION

Ref No	Qty/Mini Pack	Qty/Case
1010102	25	500



**MORRISONS**  
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)