


MORRISONS®

INFUSION SET -VENTED

I.V Set Premium



Air vent with bacteria retentive filter with cap



Latex free Y injection port for extra medication.

Product Features :

- For gravity feed only.
- Sharp piercing spike for easy insertion to IV container.
- Transparent drip chamber for easy monitoring of flow rate and priming.
- Kink resistant medical grade PVC tubing.
- Smooth roller clamp for easy flow control.
- 15 micron fluid filter for filtration of foreign particles.
- Drop rate of 20 drops /ml.
- Standard tube length-160cm.
- Single use • Disposable.
- "Sterile, individually packed"

ORDERING INFORMATION

Ref No	Qty/Mini Pack	Qty/Case
1010101	25	500



MORRISONS
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai - 600107, India.
Ph : 91 - 44 - 26155047 / 26156047.
email@morrisonslifecare.com, sales@morrisonslifecare.com
www.morrisonslifecare.com