

MOLIGUARD PLUS[®]

0.2 MICRON FLUID FILTER WITH INFUSION SET-VENTED

MOLIGUARD PLUS[®]



Air vent with bacteria retentive filter with cap



Y-site for intermediate medication



0.2 micron IV inline filter

Product Features :

- For gravity feed only.
- Infusion set with 0.2 micron IV inline filter.
- 0.2 micron IV inline filter for Entrapped air, bacterial, fungal retention.
- Kink resistant medical grade PVC tubing.
- Smooth roller clamp for easy flow control.
- Y-site for intermediate medication.
- Latex Free • DEHP free
- "Sterile, individually packed"

ORDERING INFORMATION

Ref No	Qty/Box	Qty/Case
1010110	100	400



MORRISONS
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com