


MOLIFLOW®

I.V FLOW CONTROL REGULATOR SET



- Made of medical grade polymer.
- Highly transparent kink resistant tubing.
- Flow control regulator for precise flow control of 5-250ml/hr.
- Latex free Y-injection site for extra medication site.
- Male & Female luer end for further connections.



ORDERING INFORMATION

Ref No	Qty/Box	Qty/Case
1010107	150	600



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com