

# MOLIFLOW PLUS<sup>®</sup>

I.V FLOW CONTROL REGULATOR SET WITH INFUSION SET-VENTED



- Flow control regulator with Vented Infusion Set.
- Transparent Kink resistant PVC Tubing.
- Convenient highly transparent drip chamber of drop rate of 20 drops/ml.
- Precise flow control of 5-250ml/hr.
- Latex free Y-site for intermediate medication.
- Male & Female luer end for further connections.
- Poly pouch "Sterile packing" adhering to international specification.



## ORDERING INFORMATION

Ref No	Qty/Box	Qty/Case
1010108	80	320

 **MORRISONS**  
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)