


# LIFECARE BLOOD TRANSFUSION SET NON VENTED



## **Product Features :**

- Flow rate of 20drops/ml
- Chamber with 200 micron filter
- Bulb Latex for intermediate medication
- Tip connector for needle connection
- Universal Vented spike
- 160cm length kink resistant medical grade PVC tube



## **ORDERING INFORMATION**

Ref No	Qty/Box	Qty/Case
1010106	25	400



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)