


MICRO INFUSION SET



- Flow rate of 60drops/ml
- Bulb Latex for intermediate medication
- Tip connector for needle connection
- Universal Non-Vented spike
- 150cm length kink resistant medical grade PVC tube
- Provided with 23G Needle



ORDERING INFORMATION

Ref No	Qty/Box	Qty/Case
1010104	25	500



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com