

MORRISONS BLOOD TRANSFUSION SET PREMIUM



- Flow rate of 20drops/ml
- Chamber with 200 micron filter
- Y-site for intermediate medication
- Male luer lock for needle connection
- Universal Vented spike
- 160cm length kink resistant medical grade PVC tube
- Latex Free



ORDERING INFORMATION

Ref No	Qty/Box	Qty/Case
1010105	25	400



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com