


PRESSURE MONITORING LINE



Product Features :

- Suitable for high pressure monitoring.
- Mfg. from non-toxic & non-irritant medical grade PVC.
- Highly kink resistant tubing.
- Provided with universal male luer lock & female luer lock at either end to ensure securely fits all standard equipments.
- Optional lengths are available
- Arterial / Venous lines are also available.
- "Sterile, individually packed"

ORDERING INFORMATION

Ref No	Model	Size	Qty / Box	Qty / Case
1010122	Male - Male	50cm	300	1200
1010123	Male - Female	50cm	300	1200
1010124	Male - Male	100cm	250	1000
1010125	Male - Female	100cm	250	1000
1010126	Male - Male	150cm	250	1000
1010127	Male - Female	150cm	250	1000
1010128	Male - Male	200cm	200	800
1010129	Male - Female	200cm	200	800

