

3-WAY STOPCOCK WITH EXTENSIONS



- 3-way Stopcock with a male luer lock end.
- Male luer lock end for needle / cannula connections.
- Highly transparent kink-resistant tubing.
- Available with red and blue coded tubes for arterial and venous distinction.
- Used with any of the various infusion or transfusion devices.
- Polycarbonate clear body for easy visualization of fluid flow.
- Smooth turning handle provide an efficient on/off valve and can be used in a variety of settings.
- Poly pouch “Sterile packing” adhering to international specification.

ORDERING INFORMATION

Ref No	Size	Qty/Box	Qty/Case
1010113	5cm	300	1200
1010114	10cm	300	1200
1010115	15cm	250	1000
1010116	25cm	250	1000
1010117	50cm	200	800
1010118	75cm	200	800
1010119	100cm	150	600
1010120	150cm	150	600
1010121	200cm	150	600

