

# LAPAROSUCK<sup>®</sup> PLUS


DISPOSABLE LAPAROSCOPIC TRUMPET SUCTION-IRRIGATION WITH TUBING KIT




## Specifications:

- > Sterile & Disposable
- > Soft, Flexible kink Resistant Medical Grade Tubing
- > Flexible Colour Coded Connector
- > Ergonomic trumpet valve accommodates even the most delicate hand sizes, ensuring comfortable use, minimising fatigue.
- > Easy to feel, concave/convex colour coded buttons
- > Unique detachable suction probe connection system offers fast change over probe choice.
- > Atraumatic tip design with distal vents prevents wall occlusion


**LPTF 1011**




**LPTF 2111**



**LPT 1011**



**LPT 2111**















## Insufflation Tubing With Filter



- > Prevents the risk of cross contamination
- > Luer connector for insufflation cannula
- > 0.027 micron insufflation bacterial / viral filter
- > 99.999% filter efficiency

## Our Models

Model	Design	Description	Model	Design	Description
LPT1011-5-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-5mm) and 2.5 metre tubings.	LPTF1011-5-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-5mm), 2.5 metre tubings & Insufflation tubing with 0.027 Micron Filter.
LPT1011-5-450		Suction & Irrigation Trumpet Valve with Cannula (Length-450mm, Dia-5mm) and 2.5 metre tubings.	LPTF1011-5-450		Suction & Irrigation Trumpet Valve with Cannula (Length-450mm, Dia-5mm), 2.5 metre tubings & Insufflation tubing with 0.027 Micron Filter.
LPT1011-10-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-10mm) and 2.5 metre tubings.	LPTF1011-10-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-10mm), 2.5 metre tubings, Insufflation tubing with 0.027 Micron Filter.

Model	Design	Description	Model	Design	Description
LPT2111-5-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-5mm), 2.5 metre Suction, Insufflation tubings & 2 metre Irrigation tubing.	LPTF2111-5-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-5mm), 2.5 metre Suction tubing, 2 metre Irrigation tubing & 2.5 metre Insufflation tubing with 0.027 Micron Filter.
LPT2111-5-450		Suction & Irrigation Trumpet Valve with Cannula (Length-450mm, Dia-5mm), 2.5 metre Suction, Insufflation tubings & 2 metre Irrigation tubing.	LPTF2111-5-450		Suction & Irrigation Trumpet Valve with Cannula (Length-450mm, Dia-5mm), 2.5 metre Suction tubing, 2 metre Irrigation tubing & 2.5 metre Insufflation tubing with 0.027 Micron Filter.
LPT2111-10-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-10mm), 2.5 metre Suction, Insufflation tubings & 2 metre Irrigation tubing.	LPTF2111-10-350		Suction & Irrigation Trumpet Valve with Cannula (Length-350mm, Dia-10mm), 2.5 metre Suction tubing, 2 metre Irrigation tubing & 2.5 metre Insufflation tubing with 0.027 Micron Filter.



## Morrison's Lifecare Pvt. Ltd.

No.3, Anna Street, Padikuppam, Chennai - 600 107, India.

Ph: 91 - 44 - 26155047 / 26156047.

email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com



Dealer: