


# LAPAROSUCK™

**Sterile Disposable Laparoscopic Surgery Tubing kit**




“ Providing the best solutions in single use sterile medical devices”

Model : 2111

**Specifications:**

- Sterile & Disposable
- Soft, Flexible kink Resistant Medical Grade Tubing
- Flexible Colour Coded Connector
- Prevents the risk of cross contamination
- Luer Connector for Insufflation Cannula
- 0.027 Micron Insufflation Bacterial / Viral Filter
- Models Designed as per your Convenience
- 99.999% Filter Efficiency



Model	Design	Description
1011		Sterile Suction, Irrigation, Insufflation tubing 2.5m, Colour coded connectors & Luer connector
2111		Sterile Suction, Insufflation tubing 2.5m, Irrigation tubing 2m with Spike & flow regulator, Colour coded connectors & Luer connector
3033		Sterile Suction, Irrigation, Insufflation tubing 2.5m, Colour coded connectors, Luer connector & 0.027 micron filter
3044		Sterile Suction, Insufflation tubing 2.5m, Irrigation tubing 2m with Spike & flow regulator, Colour coded connectors, Luer connector & 0.027 micron filter
Laparosuck Plus		Sterile Suction & Irrigation Trumpet Valve with Cannulae & Tubing and 0.027 micron filter
Laparosuck-GTF		Sterile Insufflation tubing 2.5m, Colour coded connector, Luer connector for insufflation cannula, 0.027 micron filter
Laparoguard		Sterile 0.027 micron insufflation Filter Bacterial & Viral filter.



## Anaesthesia

SUCTOCATH <sup>®</sup>	Suction Catheter
SUCTOCATH - TC <sup>®</sup>	Suction Catheter with thumb control
OXYGEN MASK	
OXYGEN SET	Twin Bore Nasal Oxygen Set
GUEDAL AIRWAYS	
ET TUBE	

## Urology

TOPDRAIN PLUS <sup>®</sup>	Urine Bag 2000ml with Top outlet
BOTDRAIN PLUS <sup>®</sup>	Urine Bag with Bottom outlet
LEGUBAG <sup>®</sup>	Leg Urine Bag
MOLIUMETER <sup>®</sup>	Urine Bag with Meter
MOLICATH <sup>®</sup>	Latex / Silicone Foley Balloon Catheter
KINGCATH <sup>®</sup>	Urethral Nelaton Catheter
QUEENCATH <sup>®</sup>	Urethral Catheter
ENDOFLOW <sup>®</sup>	Irrigation Set for Endoscopic T.U.R
ENDOFLOW PLUS <sup>®</sup>	Endoscopic Irrigation Set (TUR) with Outlet Drainage Connection
MOLICOND <sup>®</sup>	Male External Catheter
PAEDIATRIC URINE COLLECTING BAG	

## Gastro - Enterology

RYLES TUBE	
INFANT FEEDING TUBE	
NUTRIFEED DRIP SET <sup>®</sup>	Enteral Feeding Set
MOLIEX <sup>®</sup>	Mucus Extractor

Dealer :

## Cardio / General Surgery

CHESTOCATH <sup>®</sup>	Thoracic / Chest Drainage Catheter
MOLIVAC <sup>®</sup>	Closed Wound Suction Set
MOLISUCK <sup>®</sup>	Yankauer Suction Set
MOLI ADK <sup>®</sup>	Abdominal Drainage Kit
MOLIDRAIN <sup>®</sup>	Underwater Sealed Drainage System
CORRUGATED DRAINAGE SHEET	
LAPAROSUCK <sup>™</sup>	Sterile Disposable Laparoscopic Surgery Tubing Kit
LAPAROSUCK PLUS <sup>™</sup>	Sterile Suction & Irrigation Trumpet Valve with Cannulae & Tubing and 0.027 micron Filter
LAPAROSUCK-GTF <sup>™</sup>	Sterile Insufflation Tubing with 0.027 micron Filter
LAPAROGUARD <sup>™</sup>	Laparoscopic Insufflation 0.027 micron Filter

## Infusion / Transfusion

INFUSION SET	
BLOOD ADMINISTRATION SET	
MICRO INFUSION SET	
PRESSURE MONITORING LINES	
THREE WAY STOPCOCK	
THREE WAY STOP COCK WITH EXTENSION	
EXTENSION LINES	
MEASURED VOLUME FLUID ADMINISTRATION SET	
MOLIFLOW PLUS <sup>®</sup>	IV Flow Control Regulator with Infusion Set
MOLIGUARD PLUS <sup>®</sup>	0.2 Micron Fluid Filter with Infusion Set
PERITONEAL DIALYSIS TRANSFUSION SET	

## Miscellaneous Products

MOLICLAMP <sup>®</sup>	Umbilical Cord Clamp
NEBULIZER KIT	Mask with Nebulizer Chamber
SURGEON GOWN	
FACE MASK	
SURGEON CAP	
MOLISAFE BED	Sterile Bed Protection Sheet
MOLISAFE ADULT DIAPER	
DRAPES	
MEDICAL GRADE TUBE DIFFERENT SIZES	
(SILICONE / PVC / LATEX )	
MOLISAFE LATEX EXAMINATION GLOVES <sup>®</sup>	



CE ISO G.M.P

## MORRISONS LIFECARE PVT. LTD.

No:3 Anna Street, Padikuppam, Chennai - 600 107, India.  
Ph: 91 - 44 - 26155047 / 26156047, Telefax: 91 - 44 - 26154047  
email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com