


LAPDRAIN[®] Wound Drainage System



DESCRIPTION

- A drain bag made of superior-quality, bio-compatible siliconized PVC material.
- A drain bag provided with an efficient non-return valve prevents the reflex of fluid.
- A catheter made of kink-resistant medical-grade PVC enhances efficient drainage.
- An ergonomically designed hook with twine ensures secure hanging.
- A radio-opaque line facilitates easy and accurate placement and X-ray detection.
- Provided with an adequate catheter length of 60cm.
- Available in 10Fr, 14Fr, and 18Fr.
- Sterile, Individually packed.



ORDERING INFORMATION

Ref No.	Qty / Cover
6010165	25



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com