

Disposable FOG FREE FACE MASK – 4 LAYER



Product Features :

- 4 Layer protection with Melt-blown filter construction For Bacterial Filtration
 - BFE>99%, PFE>99% & 160mmHg Fluid resistant.
 - Soft Foam band under Nose clip for additional user comfort & prevents fogging.
 - Black fabric strips over the Nose clip deflects light glare.
 - Single-use, Disposable & Latex-free
 - Available colors: Blue & Green
- Available in : Tie-on & Earloop models.

ORDERING INFORMATION

Ref No	Model	Qty / Box	Qty / Case
9010139	TIE	100	2000
9010140	EAR LOOP	100	2000



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com