


INTRA ABDOMINAL PRESSURE MONITORING SET

(For FOLEY CATHETER)



DESCRIPTION

- Single-use for convenience and safety.
- A soft, flexible connector provides a safe grip connection with all urine drainage bags.
- 'Y' Type self-sealing, needleless connectors provide easy connection with all male luer syringes for pressure administration and sampling.
- A clear, flexible, kink-resistant, large-bore drain tube ensures uninterrupted drainage.
- Highly efficient pinch clamp-type flow regulators provide regulation of drain fluid.
- The proximal end of the tapered male connector ensures a secure connection with Foley's catheter.
- Sterile, individually packed.

ORDERING INFORMATION

Ref No.	Qty / Mini Box	Qty / Master Carton
8010120	75	300



 **MORRISONS**
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com