


MUCUS EXTRACTOR



Product Features :

- Used for aspiration of secretion from oropharynx in newborns.
- Transparent chamber for visual examination of the content.
- The capacity of the collection chamber: 25ML
- Also suitable for obtaining mucus specimens for microbiological examination.
- Smooth and round open tip for atraumatic intubation.
- Proximal end with universal connector for easy connection to the suction equipment.
- Sterile, individually packed

ORDERING INFORMATION

Ref No	Qty / Box	Qty / Case
7010118	25	500



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonlifecare.com, sales@morrisonlifecare.com

www.morrisonlifecare.com