

# MULTI LUMEN CATHETER



## Product Features :

- Manufactured from specially formulated and biocompatible PU material gives better flexibility and anti kinking performance.
- Soft and tapered tip facilitates atraumatic insertion.
- Radio opaque catheter facilitates visualization under X-ray.
- Specially designed Guide wire with soft J-tip prevents the perforation and provides good torque to ensure firm insertion.

Available Lumens : Double & Triple Lumens (Straight / Curved)

| Model                      | Size          |
|----------------------------|---------------|
| Haemodialysis Catheter Kit | 9FR - 12cm    |
| Haemodialysis Catheter Kit | 9FR - 15cm    |
| Haemodialysis Catheter Kit | 11.5FR - 12cm |
| Haemodialysis Catheter Kit | 11.5FR - 13cm |
| Haemodialysis Catheter Kit | 11.5FR - 15cm |
| Haemodialysis Catheter Kit | 11.5FR - 20cm |
| Haemodialysis Catheter Kit | 11.5FR - 24cm |
| Haemodialysis Catheter Kit | 12FR - 13cm   |
| Haemodialysis Catheter Kit | 12FR - 16cm   |



**MORRISONS**  
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)