


NASAL OXYGEN SET

Twin Bore Nasal Cannula



Product Features :

- Manufactured from soft, light weight non toxic PVC.
- Twin prong/nasal tips to ensure equal volume of oxygen to both the air passages.
- Soft funnel shaped connector facilitates easy connection to the oxygen source.
- Multichannel tube is provided to ensure continuous flow of oxygen.

ORDERING INFORMATION

Ref No	Model	Qty / Box	Qty / Case
3010136	Adult	100	400
3010137	Pediatric	100	400
3010138	Neonatal	100	400



MORRISONS
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com