

# OXYGEN FACE MASK



## Product Features :

- Moulded face mask manufactured from Non-toxic ,non irritant medical grade PVC.
- Provided with adjustable elastic strap for proper positioning of mask.
- Swivel connector for convenience of attachment to the oxygen tube.
- Multichannel tube is provided to ensure continuous flow of oxygen.
- Proximal end of tube is fitted with soft funnel shape connector for easy connection to oxygen source.

## ORDERING INFORMATION

Ref No	Model	Qty / Box	Qty / Case
3010133	Adult	50	200
3010134	Pediatric	50	200
3010135	Neonatal	50	200

