

# MOLIDRAIN-ADK

## ABDOMINAL DRAINAGE KIT



- Used for postoperative abdominal drainage.
- Collection bag of 2000ml capacity with graduated markings
- Smooth and around the distal end with large atraumatic eyes for efficient drainage.
- Radio-opaque line provided on the catheter for x-ray visualization.
- Catheter is marked at every 2cm from the last eye to ascertain the depth of placement.
- Proximal end of the tube is provided with a universal connector for easy connection to the drainage bag.
- Sterile, individually packed

### ORDERING INFORMATION

Ref No	Size	Qty / Box	Qty / Case
6010138	20fr	30	120
6010139	24fr	30	120
6010140	28fr	30	120
6010141	32fr	30	120
6010142	36fr	30	120



**MORRISONS**  
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)