

MOLIVAC[®] Closed Wound Suction Drainage Set



DESCRIPTION

- Suitable for close wound drainage under negative pressure postoperatively.
- With the option to operate one or two catheters simultaneously
- Graduated Bellow Chamber
- Non-toxic medical-grade PVC Redon drain catheters are provided with satin-smooth eyes.
- A radio-opaque line through the entire length of the catheter for radiological determination,
- Provided with a suitable-size curved knife-point stainless steel trocar.
- Sizes of the Bellow Chamber: 400 and 800 ml
- Available sizes (FR): 8, 10, 12, 14, 16, and 18 FR
- Sterile, individually packed

ORDERING INFORMATION

Ref No.	Size	Box	Master Carton
6010131	6FR	1No.	30
6010132	8FR	1No.	30
6010133	10FR	1No.	30
6010134	12FR	1No.	30
6010135	14FR	1No.	30
6010136	16FR	1No.	30
6010137	18FR	1No.	30



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com