

MOLIDRAIN-NG[®] Naso Gastric Fluid Collection Bag



DESCRIPTION

- Super smooth, kink-resistant tubing ensures unrestricted flow.
- Flexible tubing allows for patient mobility during drainage.
- Easy visualisation for when the drainage bag requires emptying.
- Graduated reading scale for volume measurement.
- Anti-reflux valve in order to reduce the contamination and infection risk and the loss of fluid in the case of capsizing.
- Provided with an outlet port for easy draining of fluid.
- A tapered step connector ensures a firm connection with the distal end of the catheter.
- A specially designed moulded handle and hook facilitate easy carrying and hanging of the bag.
- Sterile, individually packed.

ORDERING INFORMATION

Ref No.	Qty / Mini Box	Qty / Master Carton
7010119	40	160



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com