


CHESTOCATH

THORACIC DRAINAGE CATHETER WITH TROCAR



Product Features :

- Suitable for effective drainage after cardio-thoracic & thoracic surgery.
- Atraumatic & rounded open distal end with large smooth eyes for efficient drainage.
- Proximal end is fitted with a tapered tongue connector.
- Catheters are marked at every 2 cm from the last eye.
- Cross side eyes to prevent tissue aspiration.
- Radio-Opaque line and X-Ray detectable.
- Provided with sharp trocar for easy insertion
- Sterile, individually packed

ORDERING INFORMATION

Ref No	Size	Qty / Box	Qty / Case	Ref No	Size	Qty / Box	Qty / Case
6010123	10FR	20	260	6010127	24FR	20	260
6010124	12FR	20	260	6010128	28FR	20	260
6010125	16FR	20	260	6010129	32FR	20	260
6010126	20FR	20	260	6010130	36FR	20	260

