

CHESTOCATH®

THORACIC DRAINAGE CATHETER



- Suitable for effective drainage after cardio-thoracic & thoracic surgery.
- Atraumatic & rounded open distal end with large smooth eyes for efficient drainage.
- Proximal end is fitted with a tapered tongue connector.
- Catheters are marked at every 2 cm from the last eye.
- Cross side eyes to prevent tissue aspiration.
- Radio-Opaque line and X-Ray detectable.
- Provided with sharp trocar for easy insertion
- Sterile, individually packed

ORDERING INFORMATION

Ref No	Size	Qty / Box	Qty / Case	Ref No	Size	Qty / Box	Qty / Case
6010108	12 FG	20	260	6010113	28 FG	20	260
6010109	16 FG	20	260	6010114	32 FG	20	260
6010110	20 FG	20	260	6010115	34 FG	20	260
6010111	22 FG	20	260	6010116	36 FG	20	260
6010112	24 FG	20	260				

