

MOLIDRAIN[®]

Underwater Seal Drainage System



Product Features

- ♦ Underwater seal drainage system for the collection of drainage fluid from the thoracic cavity.
- ♦ Collection bag with a 1000 ml capacity and graduated markings.
- ♦ Clearly marked initial water level to ensure an underwater seal.
- ♦ Sterile, individually packed

ORDERING INFORMATION

Ref No	Qty / Box	Qty / Case
8010112-1	40	160



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com