

# MOLISUCK<sup>®</sup> Yankauer Suction Set



## DESCRIPTION

- Moulded Yankauer suction handle set
- Yankauer suction set is for convenient removal of secretion.
- Soft and kink-resistant ribbed PVC tubing of 2.5M long.
- Provided with standard connectors at both ends can connect with suction apparatus and Yankauer handle
- Yankauer suction handle is open end with 4 lateral eyes.
- Set is provided with Standard or Crown tip cannula
- Available in various lengths.
- Sterile, individually packed

## ORDERING INFORMATION

| Ref No. | Qty / Mini Box | Qty / Master Carton |
|---------|----------------|---------------------|
| 8010120 | 75             | 300                 |



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)