

EXTERNAL MALE CATHETER



Product Description

1. Morrisons silicone Male External Catheter is made of 100% silicone material.
2. The soft silicone catheter is designed to be non-irritating, non-reactive and for soft and gentle feel.
3. Clear material enables visualization of drainage flow.
4. Bio-compatible material reduces irritation for patient's comfort and safety.
5. Ergonomically shaped design enables, easily rolls over the penis.
6. The catheter is highly safe and free from any kind of allergic reactions.
7. Its proximal end is designed for safe and easy connection to all standard urine bag/ leg bag connection.
8. Individually packed ensures ready to use.

