


# ENDOFLOW<sup>®</sup>

## TRANS URETHRAL RESECTION SET



### **Product Features :**

- Y-shaped set for Endoscopic irrigation during Trans Urethral Resection of prostate gland
- Manufactured from non-toxic, non-irritant medical grade PVC
- Thumb operated clamps, helps quick and smooth change of bottles
- Proximal end is fitted with flexible latex tubing for easy connection to Endoscope

### **ORDERING INFORMATION**

Ref No	Qty/Box	Qty/Case
5010121	40	160
5010122	40	160



**MORRISONS**  
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)