

LEG U BAG[®] Leg Urine Collecting Bag



DESCRIPTION

- Manufactured from medical-grade PVC
- Super-smooth, kink-resistant, large lumen tube
- A perfect non-Return Valve to prevent backflow is provided with the bottom outlet.
- Storage Capacity: 800 ml with graduated markings
- Provided with a soft elastic strap for secure fixation
- Sterile, individually packed

ORDERING INFORMATION

Ref No.	Qty / Mini Box	Qty / Master Carton
8010120	75	300



No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : email@morrisonslifecare.com, sales@morrisonslifecare.com

www.morrisonslifecare.com