

MOLICATH[®]

SILICONE FOLEY BALLOON CATHETER



Product Features :

- Made of 100% medical grade silicone, have good biological compatibility.
- Easy insertion with the bullet-shaped tip.
- Smooth distal end provided with burr free eyes for atraumatic intubation.
- Radio opaque line for easy visualization during scans.
- Hard valves ensures easy inflation and deflation of balloon.
- Color coded sleeves for easy identification of size.
- Sterile individually packed in peelable pack.

ORDERING INFORMATION

Ref No	Model / Size	Qty/Box	Qty/Case
4010116	Two Way 6fr	10	400
4010117	Two Way 8fr	10	400
4010118	Two Way 10fr	10	400
4010119	Two Way 12fr	10	400
4010120	Two Way 14fr	10	400
4010121	Two Way 16fr	10	400
4010122	Two Way 18fr	10	400
4010123	Two Way 20fr	10	400
4010124	Two Way 22fr	10	360
4010125	Two Way 24fr	10	360

Ref No	Model / Size	Qty/Box	Qty/Case
4010126	Three Way 16fr	10	360
4010127	Three Way 18fr	10	360
4010128	Three Way 20fr	10	360
4010129	Three Way 22fr	10	360
4010130	Three Way 24fr	10	360

