

# BOTDRAIN PLUS<sup>®</sup>

## URINE COLLECTING BAG



### Product Features :

- Manufactured from medical-grade polymer.
- Super- smooth kink-resistant large lumen tube.
- Perfect Non-Return Valve to prevent backflow provided with bottom outlet.
- Specialty-designed hook or hanger facilities for carrying, handling, and holding in an upright position.
- Storage Capacity: 1800ml with graduated markings.
- Sterile, individually packed.

### ORDERING INFORMATION

Ref No	Qty / Box	Qty / Case
8010104	30	120



**MORRISON'S**  
LIFECARE PVT. LTD.

No.3, Anna Street, Padikuppam, Chennai-600107. India.

Ph : +91 96770 88178 / 72000 90887

E-mail : [email@morrisonslifecare.com](mailto:email@morrisonslifecare.com), [sales@morrisonslifecare.com](mailto:sales@morrisonslifecare.com)

[www.morrisonslifecare.com](http://www.morrisonslifecare.com)